

# HeartSine samaritan® PAD Family

## Product Comparison



**10-YEAR  
WARRANTY**



	SAM 350P	SAM 360P	SAM 500P
Compact, portable, and lightweight (1.1 kg)	✓	✓	✓
Highest level of protection from dust and moisture ingress in the industry (IP56)	✓	✓	✓
Integrated battery and electrodes for one-change maintenance	✓	✓	✓
Four-year electrode life	✓	✓	✓
Pediatric compatible	✓	✓	✓
Semi-automatic shock delivery	✓		✓
Automatic shock delivery		✓	
Audible and visual indicators for ease of use in a noisy environment	✓	✓	✓
CPR coaching with metronome	✓	✓	✓
Integrated CPR Advisor ICG-based CPR feedback			✓
Patented low escalating energy protocol (SCOPE™)	✓	✓	✓
90-minute ECG recording time	✓	✓	✓
Automatic self test	✓	✓	✓
Field upgradeable to latest AHA/ERC guidelines	✓	✓	✓
Carry case included	✓	✓	✓

### U.S./Americas

HeartSine Technologies, Inc.  
 121 Friends Lane, Suite 400  
 Newtown, PA 18940  
 Toll Free: (866) 478 7463  
 Tel: +1 215 860 8100  
 Fax: +1 215 860 8192  
 info@heartsine.com

### Europe/Rest of the World

HeartSine Technologies, Inc.  
 203 Airport Road West  
 Belfast, Northern Ireland  
 BT3 9ED  
 Tel: +44 28 90 93 94 00  
 Fax: +44 28 90 93 94 01  
 info@heartsine.com



The products described in this brochure all meet the applicable European Medical Directive requirements.

© 2014 HeartSine Technologies, Inc. All rights reserved.  
 H009-036-202-0

**Lifesaving, Pure and Simple**

[www.heartsine.com](http://www.heartsine.com)