

# NEW AED TRAINER



  
**Medico**  
*Brighter Health, Brighter Life!*

*Powerful Features  
at very affordable Price!*

# AED Trainer

## Features and functions:

- It is designed with a replaceable plug-in card which has been uploaded the languages and scenarios. To insert the card, the selected language and scenarios are installed. Once the Guidelines are changed, customers only need a new plug-in card to update the scenarios.
- What's more, with different cards, one unit can come with different languages.
- A metronome is added to the flashing lights when the unit is counting out the 2 minutes for CPR. The sound that prompts 100 compressions per minute helps trainers to properly perform CPR.
- It is completely compliant with AHA new guidelines, and pre-configured with 10 training scenarios that simulate realistic sudden cardiac arrest episodes.
- It includes the DC and AUDIO interfaces, the former can connect the adapter (that's to say, our AED Trainer can be operated by not only batteries but also the adapter), while the later AUDIO OUT can connect the speaker, it's very helpful for instructors during teaching etc.
- It is bilingual; also has Remote Control function.
- Five pairs' pads packed. Adult and paediatric pads for choice.
- There are FULLY and Semi auto scenarios selectable.



Currently there are 14 languages, if another languages are needed, that will be available.

Voice language
English (British English & American English)
French (French accent & Canadian accent)
Spanish (Spanish accent & American accent)
German
Finnish

Polish
Greek
Portuguese
Italian
Croatian
Dutch
Hungarian
Hebrew
Cantonese