

Medico Portable ECG Monitor Model: EF1800

Product Features

- Small size, Light Weight & Portable
- Standard Lead II Measurement
- 2.8 Inch TFT Color Display
- Support up to 4 different Measurement Methods
- Completes a ECG Measurement in 30 seconds
- Results include ECG Waveform, Heart Rates and Diagnosis
- Internal Memory can store up to 24 Measurements
- Results can be stored in Device or 1G SD Memory Card
- Intelligent Display Menu with User-friendly Button
- 2 Years Warranty Against Manufacturing Defects

Technical Data

- Heart Rate Range: 30 ~ 200bpm
- Heart Rate Measurement Error: $\leq \pm 5\%$
- Electrical Safety: Internally Powered, BF-Type Device
- Power Supply: 2 AAA Batteries
- Operating Temperature: 5°C ~ 40°C; Operating Humidity: <80%
- Storage Temperature: -20°C ~ 60°C; Storage Humidity: 10% ~ 90%
- Weight: Approximately 100g
- Dimensions: 125mm x 68mm x 20mm
- Included Accessories: 2xAAA batteries, 1G SD Memory Card, User Manual and Software
- Optional Accessories: ECG Lead Wires EC1000 and Disposable ECG Electrodes



For more information, please contact us at:

HTM Medico Pte Ltd

625 Aljunied Road #04-05A
Aljunied Industrial Complex,
Singapore 389836

Tel: (65) 6744-5911 Fax: (65) 6748-2718
Email: sales@htmmedico.com.sg
Website: www.htmmedico.com.sg

Insure
your whole
family's happiness
& health today!


Medico
Brighter Health, Brighter Life!



Heart disease is sudden & transient Capture it when it happens!

Concerned about palpitations, pain or shortness of breath felt by your family members or even yourself? Or that someone you value has a history of cardiovascular disease?

Now you can better monitor the condition of the heart that suffers unusual discomfort, and perhaps not fall victim to sudden cardiac death (SCD) that kills more than 1 Singaporean below 60 years old every day.

Medico

Brighter Health, Brighter Life!

Medico EF1800 is a high-tech friend that can record and store the ECG (Electro Cardio Graph) display of your heart. It is based on the hospital ECG recorder, with the added features of portability, ease of use, and timeliness that can save your heart and life.

◆ Timely Detection:

Sudden Cardiac Death is caused by an abnormal heart rhythm due to abnormal 'electrical circuitry' of the heart: it comes and goes, so that the victim is unable to record the ECG at the time when his heart is at abnormal stage. By the time the victim gets to the clinic or hospital, the root cause might be undetectable. With Medico ECG monitor EF1800, you can now timely detect and record the abnormality of your heart when you are experiencing mild chest discomfort.

◆ Diagnosis Assistance:

You can make an ECG recording with the Medico EF1800 anytime, anyplace within a minute, and the screen will display the Analysis Result, Waveform and Heart Rate. The results are accurate, reliable and can be stored, thus providing your doctor with useful information regarding your heart's condition.

◆ Prevention Before It's Too Late:

You can routinely monitor your ECG and know your heart's health for early prevention.

◆ Simple Operation:

Comes with 2.8 inch TFT Color Display and an intelligent display menu with user-friendly button.



Start Monitoring Your Heart's Health Condition
Your vital step towards a
brighter health, and a brighter life!