

PRODUCT CATALOGUE

Contents


1	Sudden Cardiac Arrest (SCA)
2	Why HeartSine?
3	HeartSafe AED in Singapore
4	HeartSine AED Models
5	HeartSine AED Accessories
6	HeartSine AED Trainer
7	CPR/AED Training Accessories
8-9	Brayden Pro Manikin: Professional Manikin with Reporting Function
10	Brayden Manikin: The First Manikin That Helps You Visualise the Effects of CPR
11	CarbonCool®: Intensive Cooling Treatment System
12-13	Cederroth: First Aid Supplies
14-15	Cederroth: Eye Wash Solutions
16	Diagnostic Tools & Personal Protection Equipments
17	Promotional Campaign

Are You Prepared for Sudden Cardiac Arrest (SCA)?




What is Sudden Cardiac Arrest?

Normal heart activity
Note the regular pattern with each beat of the heart:




Sudden Cardiac Arrest (SCA)
The heart enters a chaotic, typically ventricular fibrillation (VF) rhythm:




Problem:
The heart cannot pump blood effectively, causing the victim to collapse.


Symptoms:



Unconscious



Not breathing spontaneously



No pulse


How serious is SCA?

- Worldwide 1 million people are impacted annually
- SCA kills more people annually than lung cancer, breast cancer and AIDS combined.
- Can happen to anyone, anytime and anywhere
- Only 5-10% currently survive
- 84% of SCA events occur outside of a healthcare setting

Survival rates up to 75% have been achieved when effective CPR is applied and Automated External Defibrillators (AEDs) are readily available

How Critical is Response Time to Survival?

A Race Against Time



AEDs can increase the survival rate for SCA up to 75 percent by delivering a life-saving shock within the first few minutes of an attack. By making the decision to place AED within your company, organisation, club or school community, you are helping to save lives.

Why HeartSine?

1. **Technology Advantage**

The HeartSine AED utilizes proprietary electrode technology, advanced and stable firmware, and proprietary SCOPE™ Biphasic Technology (an escalating and low-energy waveform that automatically adjusts for patient impedance differences) to assess rhythm and recommend defibrillation if necessary.
2. **Most Portable**

HeartSine offers the smallest and lightest AED on the market, so it's the most convenient unit to use and store.
3. **Most Durable**

HeartSine AED offers the highest durability rating of any AED (including an Ingress Protection (IP) rating of IP56 and US military one-meter drop test).
4. **Most Reliable**

The HeartSine AED's 8-year warranty is the longest warranty in the industry. That includes the advanced system self-diagnostic protocols. Even HeartSine's one click PAD-Pak*—with integrated battery and electrode in a single cartridge—has one of the industry's longest electrode shelf lives.
5. **Lowest Maintenance Cost**

Unlike competitive products that require separate consumable expiration dates for customers to track and replace, HeartSine AEDs are simpler and less costly to maintain, thanks in large part to its single-unit battery electrode replacement cartridge.
6. **Lowest Cost of Ownership**

Lowest labor cost, lowest maintenance cost, field upgradable and low acquisition cost.
7. **Easiest to Use**

Designed for use by rescuers of all skill levels, HeartSine AEDs are completely self-contained units with no lid to open, complex displays or controls. Our simple two-button operation and voice/visual prompts with clear instructions guide the user through every step, including pad application, shock delivery and CPR.
8. **CPR Advisor™ Revolutionary ICG Technology.**

HeartSine® samaritan® PAD 500P provides real-time visual and audible feedback to the rescuer on the quality of compressions during a Sudden Cardiac Arrest (SCA).

Defibrillation is Easy, CPR is the Hard Part

The HeartSine® samaritan® PAD 500P provides all the advantages of a basic AED with the addition of an innovative CPR advisor. Regardless of whether it is a lay person or a professional rescuer, the AHA/ERC guidelines place a greater emphasis on the importance of CPR and advocate the prompt initiation of effective bystander CPR to significantly reduce mortality due to out-of-hospital Cardiac Arrest (CA). HeartSine can now provide a device to support these guidelines by the introduction of the HeartSine® samaritan® PAD 500P with CPR Advisor.



Easy-to-follow Visual and Voice Prompts:

1. **Built-in confidence.**

The samaritan® PAD 500P uses an Impedance Cardiogram (ICG) to access the effectiveness of the CPR.

2. **User-friendly.**

Easy-to-understand visual and voice prompts guide a user through the entire process.

Push harder/ begin CPR

Push harder

Good compression

Visual Indicators let the user know if effective CPR is being administered.

Push harder

Push Faster

Push Slower

Good Compressions

Voice prompts let the user know if chest compressions are being given correctly. Aural "clicks" help the user keep time.

HeartSine AED in Singapore

SCHOOL
Awarded the largest AED Contract to supply all government schools in Singapore. To install >2000 units of HeartSine SAM 500Ps in Primary, Secondary and Tertiary schools in Singapore.

NUS, Ngee Ann, Singapore, Republic & Temasek Polytechnics, ITE, Singapore Institute of Technology, French / Canadian / American International School, etc.

GOVERNMENT ORGANISATION
Ministry of Defense: Awarded the contract to install HeartSine AED in all Singapore Army Camps.

Singapore Civil Defense Force: Fire Bikers, Parliament of Singapore, Ministry of Health, Ministry of Manpower, Energy Market Authority, PUB, Auditor-General's Office, A*Star, Defense Science Technology Agency, etc.

HOTEL & CLUBS
Awarded period contract to supply AED to Fareast hotels, shopping centers and residences.

Marina Bay Sands, Crown Plaza Hotel, Furama Riverfront, Studio M Hotel, Orchard Hotel, Swissotel The Stamford Hotel, The Singapore Island Country Club, Aranda, Changi Beach Club, The British Club, etc.

PLACES OF INTERESTS, F&B OUTLETS
National Gallery, National Parks, Gardens by the Bay, Botanical Garden, Sports Hub, Siloso Beach Resort, Singapore Cruise Centre, Singapore Turf Club, Downtown East Resort, Aranda Country Club, Civil Service Club, Underwater World, Liang Court Shopping Center, Jumbo Group of Restaurants, etc.

PRIVATE RESIDENCE
Fareast Residences, The MarQ, Dakota Residences, Mandarin Gardens, Rio Vista, Hilltops, Martin No. 38, Queens, The Esta, The Berth by the Cove, Micasa, The Sims Ville, Varsity Park etc.

RELIGIOUS PLACES
Darul Aman Mosque, Mujahidin Mosque, Eden Bible Presbyterian Church, Heart of God, Sembawang Baptist Church, New Life Christian Church, Pelalai Buddhist Temple, etc.

BANK & TRANSPORTATION
DBS Bank, OCBC Bank, May Bank, Julius Bar, BSI, UBS, DHL, FedEx, SMRT Bus Depot, SMRT Taxi, etc.

CONSTRUCTION, MARINE, OFFICE, FACTORY, CLINICS, ETC.
Kajima Construction, Sato Kohyo, Hallin Marine, Schlumberger, Canyon, Eastman Chemical, Roche, TUV, Novartis, Eastman, Natsteel, JTC Corporation, Canyon Offshore, Nestle, Pfizer, Avago Technologies, IBM, PSB, Aesthetic & Reconstructive Centre, Tzu Chi Free Clinic, etc.



HeartSine AED Models

Product Comparison



	SAM 350P	SAM 360P	SAM 500P
Compact, portable and lightweight (1.1kg)	✓	✓	✓
Highest Level of protection from dust and water ingress in the industry (IP56)	✓	✓	✓
Integrated battery and electrodes for one-change maintenance	✓	✓	✓
Four year battery/electrode life	✓	✓	✓
Pediatric compatible	✓	✓	✓
Semi-automatic shock delivery	✓		✓
Automatic shock delivery		✓	
Audible and visual indicators for ease of use in a noisy environment	✓	✓	✓
CPR coaching with metronome	✓	✓	✓
Integrated CPR Advisor ICG-based CPR feedback on Force and Rate (Push harder, Push faster, Push slower, Good compression)			✓
Proprietary low escalating energy protocol (SCOPE™)	✓	✓	✓
90-minute ECG recording time	✓	✓	✓
Automatic self test	✓	✓	✓
Field upgradable software to latest AHA/ERC/NRC guidelines	✓	✓	✓
Carry case included	✓	✓	✓
8 Years Warranty	✓	✓	✓

AED Accessories

AED Cabinets & Cases



Floor Standing Cabinet Table/Wall mount Cabinet



Rigid Watertight Case




AED Rescue Backpack

AED Signage



2-Way AED Sign 90° Double-sided AED Sign Flat Wall AED Sign



AED Directional Sign

Battery & Defibrillation Pads



Pediatric Cartridge (4 Year Shelf Life)



Adult Cartridge (4 Year Shelf Life)

Other Accessories



1st Aid Respond Kit Consist of pouch, CPR mask, shaver, scissor, glove, towel



USB Data Cable

Picture for reference only. Subject to change without prior notice.

HeartSine AED Trainer

Provides Realistic Training Demonstration and Superior Retention

The easy-to-use samaritan PAD Trainer with remote control guides users through simulated analysis, energy delivery and prompted Cardio Pulmonary Resuscitation (CPR) intervals with the look and feel of a live samaritan Automated External Defibrillator (AED) without the actual charge and discharge of an electric shock.

The HeartSine® samaritan® AED Trainer package includes remote control with AAA batteries, Trainer-Pak, battery charger, soft carry case with two zipper compartments, user manual, and three pairs of training electrode pads.



Features:

- **Cable rewind mechanism.** Easy replacement of pads and cables after each session reduces time spend resetting system.
- **Rechargeable Battery.** Operates up to seven hours-a full day of training-on a single charge. Recharges in just four hours.
- **Battery status indicator.** LEDs indicate charging status (while not in use) or one hour of remaining battery life (While in use).
- **Pause Button.** Temporarily stops, then resumes, the training scenario to enable the trainer to offer additional explanation and support to the class.
- **Volume Configurable.** Four selectable volume levels with default to last selection.
- **Selectable Training Scenarios.** Six pre-programmed AHA/ERC scenarios.
- **Training Pads.** Three pairs of reusable, replaceable electrode pads included. Compatible with all training manikins,
- the new electrode pads are now easier to place back on the pad liner.
- **Upgradeable.** Updates to Trainer software available on the HeartSine website.
- **Language Options.** Twenty-nine language options enable the Trainer to provide instructions in the local language.
- **Warranty.** 2 Years against manufacturing defects.

Simulate the Full Range of "Real Life" Scenarios

The multi-function remote control enables you to select among six AHA/ERC scenarios for live demonstrations or training. One remote control can even operate multiple trainer units at a time (within a range of 3-5 metres/10-15 feet). The Trainer is also compatible with any CPR manikin system.

Realistic

Voice and visual prompts and CPR Advisor match those of the HeartSine® samaritan® AED defibrillator. Likewise, size and shape are the same. Unique labeling and membrane color prevent inadvertent use in an actual emergency.

Appropriate for All Skill Levels

The trainer can be used by businesses to train employee first responders on the use of the HeartSine® samaritan® AED. In addition, it meets the needs of professional trainers, either teaching first-time lay users or refreshing experienced emergency responders.

Cost-Effective

The Trainer is powered by a rechargeable battery. Electrode pads are cost-effective and reusable, and mimic the PAD-Pak™ battery/pad cartridge unique to the HeartSine® samaritan® AED defibrillator.

Medico AED Trainer

Cost Efficient & Multi-Functional AED Trainer


This Medico AED trainer has all the functions needed to teach your students about AED operation – without breaking your budget. This lightweight, affordable trainer offers bilingual scenarios (English and Chinese), a remote control, 3 sets of adult pads, 2 sets of child pads, connectors, all packed in a nylon carrying case. It runs on AA batteries or the AC adapter and has an audio port for connecting speakers for bigger crowd.



CPR / AED Training Accessories

HeartSine AED Trainer Accessories

- Trainer Cartridge
- Training Defibrillation PAD



Medico AED Trainer Accessories

- Training Defibrillation PAD



Training Face Shield



CPR/AED Training Manikin

- Brayden CPR Manikin
- Prestan CPR Manikin



AED Product Refresher Training & Certified CPR/AED/First Aid Courses Available. Call us for details.

ACLS/BCLS Professional Manikin with Reporting Function

The Brayden Pro provides healthcare professionals with quantitative data driven feedback. This can be used to enable students to improve their CPR performance when a data driven need is required. Brayden Pro provides students with real time feedback during CPR training and also enables detailed post CPR training debriefing by providing quantitative data of key CPR parameters as outlined in the ERC 2015 Guidelines such as compression depth, compression rate, hand position and compression release. Feedback on hands off time and ventilation performance is also provided from the Brayden Pro manikin.

"The Bluetooth link between the impressive Brayden Pro manikin and the app on an android tablet greatly enhances the training value of the system. It gives not only an immediate guide to the quality of compressions, but also offers trends and a permanent record of what is being achieved. It would be hard to rival as a training system."

PROFESSOR DOUGLAS CHAMBERLAIN



BRAYDEN^{PRO}

"The Brayden Pro is a breath of fresh air in the world of resuscitation. The pedagogical feedback mechanism through direct visual LED lights & the link to the Android tablet platform gives a superb development in auditing the student and giving them a clear understanding of direct and future learning. The equipment is easy to use and has not been overcomplicated by unnecessary bells and whistles."

BOB FELLOWS,
PRE-HOSPITAL HEALTHCARE EXPERT




ADVANCED CPR COACHING

- Quality Score using algorithm determined by leading European Opinion Leaders
- Can be configured for AHA & ERC (2015) guidelines
- Can be set for differing responder skill level - Healthcare Professional or Layperson
- Can be set for one or two responder scenario
- Can be set for ventilation only and compression only
- Tested & Evaluated by European in-hospital & pre-hospital resuscitation experts.

"The Brayden Pro is a cost effective CPR training manikin for healthcare professionals. Clinical staff require more precise skills training. This upgraded Brayden Manikin with tablet can deliver accurate real time and post event feedback. The tablet software displays detailed information of both compressions and ventilation in a very user friendly way. This allows precise individualised feedback. This is a welcome addition to the developing Brayden Manikin."

ROB MORRISON,
IN-HOSPITAL PROFESSIONAL



REAL TIME CPR FEEDBACK

- Up to 6 manikins can be monitored at one time.
- Connected via bluetooth.
- Intuitive app & feedback screens.

OBJECTIVE ANALYSIS & DEBRIEFING

- Timeline of assessed performance measured through 8 key CPR metrics.
- Quantative analysis of key CPR metrics for debriefing.
- Can be configured to AHA and ERC guidelines.

OVERALL CPR QUALITY SCORE CALCULATED & PRESENTED OVER TIME

- Objective CPR metrics available as complete event & by individual CPR cycle.
- Android based- data can be exported.


The First Manikin That Helps You Visualise the Effects of CPR

BRAYDEN

CPR QUALITY INDICATING LAMPS

These lamps are only light when the compression depth and the speed coincide correctly with the given default value. When either one of the two i.e. the compression depth or the speed is not correct then the lamps will fail to light.

The state of the art technology integrated into "BRAYDEN" makes it possible to train all levels of trainee, making it simple and easy to learn the correct way to perform CPR effectively. The trainee can learn simply by focusing attentively on the indicating lamps and how they respond to their CPR performance.



BLOOD CIRCULATION INDICATING LAMPS

These lamps provide visualisation of the blood flow from the heart to the brain which trainees can easily observe while performing CPR. The lamps will illuminate in direct proportion to the depth and speed of the trainees compressions providing feedback throughout their CPR practice sessions.

COMPRESSION INDICATING LAMP

These lamps directly relate to depth and speed of compressions. In other words, these compression indicating lamps across the chest light up in direct proportion with the compressions depth and speed as with the Blood Circulation Indicating Lamps.

For example, when the CPR is correctly performed all the lamps across the chest are fully lit, however the lamps will only light up in proportion with the quality of CPR.

"BRAYDEN" is the first CPR Manikin ever to display the virtual flow of blood from the heart to the brain during the CPR

✓ Builds Confidence

✓ Promotes Understanding

✓ Aids Knowledge Retention

Compression Cooling System, Dry & Portable

CarbonCool®

The CarbonCool® intensive cooling system is dry, portable, and allows for quick and easy application to the human body for quick and easy application to the human body for compression cooling. The CarbonCool® system allows a snug fit to any body type and is also suitable for multiple use.



ACTIVE & RECOVERY COOLING

Firefighters

Rapid body cooling is applicable for firefighters to reduce core body temperature during and after training and/or operation. Active cooling is proven to reduce mean core temperature, heart rate and leukocyte count of Firefighters, preventing heat stress and possible injuries.

Sportsmen

In Sports, cold and compression (in R-I-C-E principles) are routinely applied immediately after acute injury to alleviate pain, reduce swelling and speed functional recovery.



MEDICAL APPLICATION

- ✓ Reduce Inflammation
- ✓ Reduce Swelling
- ✓ Reduce Pain
- ✓ Reduce Muscle Spasm
- ✓ Reduce Cardiac Output
- ✓ Reduce Oxygen Consumption
- ✓ Improve Cardiovascular Stability (circulatory disorder)
- ✓ Improve Serum Lactate Clearance
- ✓ Targeted Temperature Management
- ✓ Continue Temperature Management beyond 24 hours
- ✓ Out-of-Hospital Cooling



COMFORT COOLING (HAZMAT)

For personnel donning protective wear such as the Hazmat suit, our cooling system is able to provide sustained comfort and heat stress reduction.

Reduces Risk of Heat Injury

Continuous cooling aids to extract excess heat and maintain a balanced temperature within the protective suit, preventing overheating.

During Operations Cooling

Cooling during operations helps keep the wearer comfortable within the suit, boosting alertness and work capacity.



Cools faster, lasts longer. Designed for cold compression and multiple usage across many verticals

First Aid Supplies

Cederroth First Aid Station

The world's smallest emergency department has now opened. The Cederroth First Aid station is more hygienic and can be adapted to deal with different risks and healthcare needs.

Cederroth First Aid Station

Contents

- 2 Cederroth 4-in-1 Bloodstoppers
- 3 Cederroth 4-in-1 mini Bloodstoppers
- 1 Triangular Bandage (2-pack)
- 2x20 Salvequick Wound Cleansers
- 1 Salvequick Dispenser including 40 Fabric and 45 Plastic adhesive bandages (sterile)
- 1 First Aid Instructions
- 1 Key
- Dimensions B 29x H56 x D12cm

REF 490940

Hygienically Clean!

A common problem at many workplaces is that First Aid equipment can get dusty and dirty. This new design includes a protective door, which keeps your equipment clean and ready to use.

Clear display

The first Aid instructions on the station are easy to follow and there are clear instructions on product packages.

No unnecessary waste

The adhesive bandage refills are locked in the dispenser to make sure they do not disappear. The wound cleanser dispenser is also locked. It is easy to remove the empty refills using a special key.

Space for extra compartments

There are two additional storage compartments, helping you to adapt your station to the specific needs of your workplace. You could use these compartments for Cederroth Eye Wash Pocket Model (REF 722110) or Cederroth Soft Foam Bandage (REF 676100) for example.

Adhesive bandages that you pull downwards

In our new dispenser, you simply pull the adhesive bandages downwards, which avoids getting any dirt or blood on any adhesive bandages.

CEDERROTH



First Aid Panel Mini

This is a small panel that is suitable for small spaces or for work vehicles such as lorries and vans.

Contents:

- 1 Salvequick Dispenser including 40 Fabric and 45 Plastic adhesive bandages (sterile)
 - 1 Cederroth 4-in-1 Bloodstopper
 - 3 Cederroth 4-in-1 mini Bloodstoppers
- Dimensions: B25 x H40 x D5 cm.

REF 191460



Plaster Dispenser

Innovative plaster dispenser allows you to pull adhesive bandages downwards and put on with only one hand. This avoids you getting any dirt or blood on any other bandages. No unnecessary wastage as refills are locked in the dispenser.

Contents:

- 40 Fabric Sterile adhesive bandages
 - 45 Plastic Sterile adhesive bandages
 - 1 key
- Dimensions: B23 x H12 x D5.5 cm

REF 490770



Salvequick Plastic Adhesive Bandages

Each refill holds 45 sterile adhesive bandages (27 pcs measuring 72 x 19mm and 18 pcs measuring 72 x 25 mm) 6 refills/box.

REF 6036 CAP



Salvequick Blue Detectable

Sterile, extra soft and flexible. Detectable by most metal detectors used in the food industry. Each refill holds 35 sterile adhesive bandages (21 pcs measuring 72x19mm and 14 pcs measuring 72 x 25 mm) 6 refills/box.

REF 6735 CAP



Cederroth Soft Foam Bandage

Cederroth Soft Foam Bandage is a tight-fitting elastic bandage, giving users maximum freedom of movement and it will stay in place in water. Contains latex, Roll 2 metres x 6 cm, pack of 2.

REF 676100

Glue-free. Self-adhesive. Water-resistance. Elastic



Cederroth 4-in-1 Bloodstopper

Sterile universal dressing with four functions: pressure pad, protective dressing, temporary support dressing or burn dressing. Contains 1 dressing (14x23cm) and 2 elastic bandages (10cmx3m). Clear instructions on the inside of the pack.

REF 1910



Salvequick Wound Cleanser

Salvequick Wound Cleansers do away with cotton wool and bottles of antiseptic. Contents: 20 Salvequick Wound Cleansers (11x12cm) with 0.9% sodium chloride solution.

REF 323700



Cederroth Triangular Bandage

This can be used to support an injured arm or to keep a dressing in place on a head injury. Contains 2 triangular bandages and 6 safety pins.

REF 1889

Cederroth Eye Wash

Using the right Equipment Improves the Chance of saving your eyes

When alkalis or acids get into the eye, the pH value of the eye changes and corrosion occurs. Alkaline substances, for example **lye, cement, washing soda and dishwasher detergent**, are the reasons behind most corrosive injuries. In the event of an accident, the first thing you need to do is to quickly start to flush away the harmful chemicals and restore the pH value of your eye. This is why Cederroth's eye wash is buffered. It works as a shock absorber and restores the normal pH value of the eye more effectively than normal sodium chloride solution.



It's a good idea to store the bottle in the wall bracket, so it has a fixed place and everyone knows where it is.

Generous flow of liquid for 1.5 minutes

The bottle is designed to rinse the eye with a generous flow of liquid for around 1.5 minutes. Make sure you have more bottles at hand to rinse for a longer period of time, at least 15 minutes.

4.5 year shelf life!

Cederroth Eye Wash has one of the longest shelf lives on the market for these kinds of products, which saves money and improves safety

Proven to be more effective for alkali and acid splashes

Cederroth's buffered liquid helps to bring the pH value back to neutral again, unlike a normal sodium chloride solution or water, which simply dilutes and rinses.



The bottle is easy to open

The bottle opens easily by firmly twisting the green dust cover either clockwise or anti-clockwise. If the bottle is in the wall bracket, it will open as soon as you twist it out of the bracket. The size of the bottle makes it easy to hold.

The integrated eye cup guides the liquid towards the eye without any unnecessary waste

The eye cup helps to keep the eye open during rinsing. A colleague can check that the eye is open by looking through the transparent cup.

Clear graphical instructions on the bottle showing you what to do


100% phosphate-free

If you have an accident with a substance that contains calcium (cement, concrete or lime), you can safely rinse with Cederroth Eye wash as it is 100% phosphate-free. Calcification may occur with eye rinses that contain phosphates — the phosphate in the liquid reacts with the calcium, forming an insoluble white deposit which can further damage the eye.

CEDERROTH

CEDDEROTH EYE WASH IS PROVEN MORE EFFECTIVE FOR ALKALI AND ACID SPLASHES.

When alkalis or acids get into the eye, the pH value of the eye changes and corrosion occurs. Cederroth's borate-buffered liquid helps to bring the pH value back to neutral again, unlike an unbuffered eye wash, which simply dilutes and rinses. The eye has a neutral pH value (7). Acids have a low pH value and alkalis have a high pH value. This is why Cederroth Eye Wash has been recommended by eye specialists in a study that compared the effect of different rinsing solutions on alkalis¹.



1) Rihawi et al., Graefes' arch Clin. Exp.ophtamol. (2006) 244:845-854

Complete Cederroth First Aid Solutions

It is particularly important to have first aid in place for handling emergency accidents in workplaces with a high risk of eye injuries, for example **construction sites, warehouses, catering kitchens, industry and laboratories**. Discover the entire range at www.firstaid.cederroth.com



Diagnostic Tools

Medico Portable ECG Monitor EF1800



Portable ECG Monitor for daily home use

- ECG & Heart rate values can be easily measured by layperson and doctors
- For home, leisure, sports and work
- Easy and simple one-hand operation
- Lightweight and compact
- Complete with standard accessories, analysis software and 2 years warranty
- CE Certified, patented



Heart disease is Sudden & Transient, capture it when it happens!

Pulse OxiMeter



A pulse oximeter is a medical device that indirectly monitors the oxygen saturation of a patient's blood (as opposed to measuring oxygen saturation directly through a blood sample) and changes in blood volume in the skin, producing a photoplethysmogram.

Thermometer



- Contact Thermometer
- Non-contact Thermometer

Personal Protection Equipments



Head Protection

- Safety Helmet
- Protective Face Shield
- Safety Goggle



Hand Protection

- Safety Gloves
- Fire Retardant Glove



Leg Protection

- Safety Boot
- Firemen Boot



Body Protection

- Fire Retardant Garment
- Safety Apparel



Communication

- Loud Hailer
- Walkie-Talkie



Fire Protection

- Fire Blanket
- Fire Extinguisher

Forward Hearts

Enables individuals who have survived a sudden cardiac arrest event as a result of a HeartSine defibrillator to save another's life by donating a samaritan® Public Access Defibrillator to the charity or organization of their choice. To celebrate the life-saving event, the donation will be made in the survivor's name, in HeartSine's name, and in the distributor's name.

HOW THE PROGRAM WORKS:

When HeartSine AED device is used to save a life, the distributor works to connect everyone involved - the company site, the saved individual, and HeartSine - to begin the donation program. The saved individual can choose to make his/her donation to the charity anonymously and confidentially, or as part of a public announcement.

Read more about our Forward Hearts program at www.heart sine.com/forwardhearts



Free Pad-Pak™ Program



If you have used a HeartSine samaritan® PAD during a Sudden Cardiac Arrest (SCA), you're entitled to a free Pad-Pak when you provide the Saver EVO file, along with a completed Customer Event Form.

FREE PAD-PAK RULES:

- Within 72 hours after the use of the AED during an SCA, contact us with the following information:
 - Download the event data from the AED
- Complete and submit the HeartSine Customer Event Form with information about the patient
- Once the information has been received and approved by HeartSine, you will be entitled to the free Pad-Pak.

It's Lifesaving. Pure and Simple.

“HTM” stands for “ HELP the MASS ”

Our Mission is to be the leading supplier of good quality and affordable Portable / Home Medical Devices and Emergency Life Saving products for lay people.

Our business Philosophy & Vision is to educate the masses on how to take charge of our own health and embrace a healthier life thus achieving our Slogan

“Brighter Health, Brighter Life!”

HTM Medico Pte Ltd

625 Aljunied Road #04-05A
Aljunied Industrial Complex Singapore 389836

 **Tel:** (65) 6744 5911  **Fax:** (65) 6748 2718

 **Email:** sales@htmmedico.com.sg  **Web:** www.htmmedico.com.sg



HeartSine®

CarbonCool®

CEDERROTH

BRAIDEN