

# FREE Pad-Pak™ program

## Sudden cardiac arrest? Tell us.

If you have used a HeartSine samaritan PAD during a sudden cardiac arrest (SCA) in a public setting, it's critical that you TELL US. To do this, provide HeartSine with the Saver EVO file that contains comprehensive event data about the event and a completed Customer Event Report form. This information is used by HeartSine to measure the success of our automated external defibrillators (AEDs) and to continually improve our products to ensure optimum results. When you provide the Saver EVO file, along with a completed Customer Event Report form for a qualifying\* event, we'll send you a free Pad-Pak for your efforts!



## Two simple steps

- 1** Immediately after a HeartSine AED is used during an SCA, download the event data file (.evo) from your AED.
- 2** Complete and submit the following to [aedevent@stryker.com](mailto:aedevent@stryker.com):
  - HeartSine Customer Event Report form (with no identifying patient information)
  - Event file (.evo)

If the event qualifies\* and the information provided is complete, a free Pad-Pak will be shipped to you. Please allow two business days for review of the event data.

The Saver EVO software and manual, and the Customer Event Report form are available at [heartsine.com](http://heartsine.com).

### \*Qualifications

- Event must be a sudden cardiac arrest to qualify. (Event is reviewed by the HeartSine Clinical team whose decision is final.)
- Please refer to [heartsine.com](http://heartsine.com) for the complete list of requirements to qualify for the Free Pad-Pak program after your device has been used in a sudden cardiac arrest.

**For further information, please contact your Stryker representative or visit our website at [strykeremergencycare.com](http://strykeremergencycare.com)**

## Emergency Care Public Access

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: HeartSine, Pad-Pak, samaritan, Saver EVO, Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.


Date of Issue: 06/2021

Made in U.K.

H009-032-211-AE EN-UK

HeartSine SAM 500P is not available for sale in the U.S.

Copyright © 2021 Stryker



HeartSine Technologies Ltd.  
207 Airport Road West  
Belfast, BT3 9ED  
Northern Ireland  
United Kingdom  
Tel +44 28 9093 9400  
Fax +44 28 9093 9401  
[heartsinesupport@stryker.com](mailto:heartsinesupport@stryker.com)  
[heartsine.com](http://heartsine.com)

Stryker European  
Operations B.V.  
Herikerbergweg 110  
1101 CM Amsterdam  
Netherlands  
Tel +31 (0)43 3620008  
Fax +31 (0)43 3632001