

# HeartSine® Samaritan PAD family



	HeartSine SAM 350P	HeartSine SAM 360P	HeartSine SAM 500P with CPR Advisor™
ICG-based real-time feedback on CPR rate and force (“Push harder,” “Push faster,” “Push slower,” “Good compressions”)			✓
CPR coaching with metronome	✓	✓	✓
Semi-automatic shock delivery	✓		✓
Fully automatic shock delivery		✓	
High level of protection from dust and water ingress (IP56)	✓	✓	✓
Integrated battery and electrode for one-change maintenance	✓	✓	✓
4-year electrode and battery life	✓	✓	✓
Low energy escalating protocol (SCOPE)	✓	✓	✓
Audible and visual metronome for CPR rate guidance	✓	✓	✓
Compact, portable, and lightweight 1.1 kg	✓	✓	✓
Pediatric compatible	✓	✓	✓
Automatic self test	✓	✓	✓
Carry case included	✓	✓	✓
Size: 20 cm x 18.4 cm x 4.8 cm			